

START HERE GET TO KNOW THE BUTTONS & ATTACHMENTS

Please carefully review all Warnings, Precautions, and Contraindications in User Guide prior to first use.



EXFOLIATION ATTACHMENT:

INTENSITY BUTTON

Short press to change suction level.

POWER BUTTON

Short press to pause and resume.




DEPUFFI ATTACHMENT:

INTENSITY BUTTON

Short press to change temperature intensity.

POWER BUTTON

Short press to switch between **INSTACHILL** and **INSTAHEAT** settings.



EXFOLIATION BENEFITS



ENHANCED GLOW

Radiant, lit-from-within glow.



SUPERCARGED EXFOLIATION

Combines AHA/BHA and physical exfoliation for brighter, more radiant skin.



POWERED DETOX

Clearer, less visible pores.



BOOSTED HYDRATION

Leaves skin hydrated for a smoother feel.



DEPUFFI BENEFITS



INSTACHILL

ADVANCED DEPUFFING

Reduces appearance of dark circles and leaves under eyes refreshed and visibly less puffy.

SCULPTED APPEARANCE

Visibly enhance the appearance of cheekbone and jawline definition.

RAPID SOOTHING

Cooling relief for instantly calmer soothed skin.



INSTAHEAT

ACTIVATED RADIANCE

Massage-driven circulation for a healthy, radiant look.

SOFTENING PREP

Primes skin for extraction by loosening surface impurities.

FACIAL DE-STRESS

For muscle tension relief.



CLINICALLY TESTED SKINCARE

Dermatologist-Tested Korean skincare formulas.
Developed for optimal results when using FacialPro Glow device.



DERM DETOX

AHA + BHA Exfoliating Gel

Gently softens and dissolves excess dead skin cells and sebum sitting on the surface of skin.

APPLY



HYDRO INFUSE

BHA Hydrator

Hyaluronic acid infuses skin with moisture.

INFUSE



Pro Tip: Don't forget to use both Wide and T-Zone tips when using your Hydro Infuse for best at-home results.

TO REPLENISH SKINCARE, VISIT [SHARKBEAUTY.CO.UK](https://sharkbeauty.co.uk)

Only use this skincare with FacialPro Glow system.



SCAN TO
REPLENISH

Shark

FACIALPRO GLOW

+ DePuffi

QUICK START GUIDE



SCAN FOR
HOW-TO VIDEOS



DAILY DEPUFFING ROUTINE

Depuff to sculpt and boost radiance daily, in just 4 minutes.

GET STARTED

- Clean your face to remove any dirt and oil.
- Connect DePuffi Attachment.
- Apply your favourite serum, face oil or moisturiser for best glide.



INSTACHILL / 4 MIN. SESSION

- Turn the device on.
- Device starts at temperature level 3. Start there and dial down to your preferred temperature setting.
- Glide over face, neck, jawline, and shoulders to reduce puffiness for a sculpted look and refreshed skin.



INSTAHEAT / 4 MIN. SESSION

- Turn the device on.
- Short press to switch into InstaHeat mode.
- Device starts at temperature level 3. Start there and dial down to your preferred temperature setting.
- Massage and glide over face, neck and jawline to boost circulation, encourage lymphatic flow and boost skin's radiance.



DEFINE AND DEPUFF

INSTACHILL

- Move in slow upward and/or outward strokes to reduce puffiness for a defined jawline and cheekbones that look sculpted.
- Use under eyes to reduce the appearance of dark circles for an instantly refreshed and awakened look.



BOOST RADIANCE

INSTAHEAT

- Massage in slow upward and/or outwards strokes to boost circulation for more radiant skin.
- Use in downward motion on the neck for lymphatic flow.



WEEKLY PRO GLOW ROUTINE

Smooth, brighten, and hydrate with a weekly skin refresh and glow in just 9 minutes.



DEPUFFI RECOMMENDATION: Prior to exfoliation, prep skin with **INSTAHEAT** to loosen surface impurities.

1

APPLY DERM DETOX / 3 MIN. SESSION AHA + BHA gently dissolves impurities

- Apply a thin layer of **Derm Detox** all over face, avoiding eyes and mouth.
- Leave on for 3 minutes and proceed to next step.



2

EXFOLIATE & EXTRACT / 3 MIN. SESSION Clears pores and lifts impurities.

- Fill tank with cool water to “MAX” fill line.
- Attach **Wide Tip** and turn device on.
- Adjust intensity and glide over skin as directed. Pause device and attach **T-Zone Tip** for targeted care. Device will resume at time you left off
- Device pauses after 3 minutes—move to next step.
- Rinse off any remaining **Derm Detox** from face with water.

FILL



CHOOSE



EXFOLIATE



3

HYDRO INFUSE / 3 MIN. SESSION Locks in moisture and refreshes skin.

- Fill tank with **Hydro Infuse** to the top of the “MAX” fill line.
- Select your **tip** and intensity level as Step 2 and glide over skin as directed.
- Once you have completed your Weekly Pro Glow Routine, apply moisturiser and/or SPF
- If using DePuffi on **INSTACHILL** as recommended below, apply moisturizer beforehand and SPF after completing the step.

FILL



CHOOSE



INFUSE



DEPUFFI RECOMMENDATION: After extraction, refresh skin by applying your moisturiser, using **INSTACHILL** to soothe for instantly calmer skin and reduced redness. Complete routine by applying SPF if it's daytime.

SELECT EITHER THE NORMAL TIPS OR THE GENTLE TIPS BASED ON YOUR PREFERENCE.

CHOOSING A SET OF TIPS

All tips were designed with performance in mind.

Dermatologist-tested to deliver clean pores and hydrated, glowing skin.

NORMAL



stronger feeling of suction
FOR MOST SKIN TYPES

CHOOSE
A SET



GENTLE



more mild feeling of suction
FOR SENSITIVE SKIN TYPES

DO NOT HOLD THE TIPS STILL ON THE SKIN. KEEP THEM MOVING.

2MIN




HOW TO GLIDE WITH **WIDE TIP**

For use on cheeks, jawline, and brow bones.

We recommend you use the Wide Tip for 2 out of the 3 mins from both Step 2 and 3

- Place the tip on your cheek, jawline, or upper forehead and press gently to create a seal.
- You'll feel suction and hear a soft "kiss" with each stroke.
- Move in short, smooth strokes about 2-5 cm long.
- Use your free hand to gently stretch skin for a smooth glide.



1MIN




HOW TO GLIDE WITH **T-ZONE TIP**

For use on T-zone: nose, chin, and around eyebrows.

We recommend you use the T-Zone tip for 1 out of the 3 mins from both Step 2 and 3

- Place the tip on your nose, chin, and between eyebrows and press gently to create a seal.
- You'll feel suction and hear a soft "kiss" with each stroke.
- Move in short, smooth strokes about 2-5 cm long.
- Use your free hand to gently stretch skin for a smooth glide.



The combination of both **Wide** and **T-Zone Tip** is what will help to achieve optimal results across all areas of the face.